


PERFUSION SYSTEM FOR MESENTERIC BEDS

Study the effect of vasoactive substances and electrical stimulation



FEATURES & BENEFITS

- » Up to 4 mesenteric beds may be studied
- » Perfusate flow rate from 5 to 20 ml/min
- » Several small cannulas to suit different artery sizes
- » Easy to manipulate stainless steel grids
- » Adjustable bath heights and volumes
- » Adjustable carbogen supply
- » Easy maintenance
- » Compact
- » Optional electrical stimulation through stainless steels beds
- » One protocol for each organ (different flow rate from one bath to another, simultaneous or staggered perfusion...)



PRINCIPLES

- » Each organ lies on a stainless steel grid and bathes in a buffer solution
- » The organ is perfused at a constant flow rate via its artery
- » Buffer flow rate 5 to 10 times greater (max 150ml/min) than perfusion flow in order for perfusion to fully infuse bath liquid
- » Pressure is measured just upstream of the cannula, reflecting changes in vascular tone
- » Buffer solution maintained at the right temperature with heated water, circulating behind the bath walls
- » Carbogen is released directly into the buffer solution



.../...

SETUP

- » 4 independent tissue baths (90 x 60 x 20 mm)
- » Flow regulator (peristaltic pump)
- » Heater/pump for heating circuit
- » Amplifier

ACQUISITION & REAL-TIME ANALYSIS

The system can be used with emka TECHNOLOGIES's iox2 software:

- » Perfusion pressure measured upstream of each mesenteric bed is acquired
- » iox' analyzer generates an average value based on user-defined intervals
- » External devices control
- » Protocols, to automate the experiment

